


Kinex

VEGA

ACL Fixation System



OSTEONIC
The Medical Explorer

BIO Interference Screw


Overview

The Bio interference screw is a composite of 70% PLGA(lactic-co-glycolic acid) and 30% micro β -TCP(Beta-tricalcium phosphate)


Feature benefits

- ✿ High torsional strength
- ✿ Blunt thread design to prevent tendon grafts from damage
- ✿ Tapered design for easy insertion


Type I

Item	Item Code	Description
	CIS0725-BC	Ø7.0*25mm
	CIS0730-BC	Ø7.0*30mm
	CIS0825-BC	Ø8.0*25mm
	CIS0830-BC	Ø8.0*30mm
	CIS0925-BC	Ø9.0*25mm
	CIS0930-BC	Ø9.0*30mm
	CIS1025-BC	Ø10*25mm
	CIS1030-BC	Ø10*30mm

Type II

Item	Item Code	Description
	CIS0823-BC	Ø8.0*23mm
	CIS0828-BC	Ø8.0*28mm
	CIS0923-BC	Ø9.0*23mm
	CIS0928-BC	Ø9.0*28mm
	CIS1023-BC	Ø10*23mm
	CIS1028-BC	Ø10*28mm

Instrument

Item	Item Code	Description
	ACL-0329-DRIVER-1	Bio interference screw driver


PEEK Interference Screw

Overview

The Peek interference screw is made of radiolucent peek material

Feature benefits

- ✦ Strong, rigid fixation
- ✦ Radiolucent material
- ✦ Enhanced durability

Item	Item Code	Description
	PIS0725-PK	Ø7.0 x 25 mm
	PIS0730-PK	Ø7.0 x 30 mm
	PIS0825-PK	Ø8.0 x 25 mm
	PIS0830-PK	Ø8.0 x 30 mm
	PIS0925-PK	Ø9.0 x 25 mm
	PIS0930-PK	Ø9.0 x 30 mm
	PIS1025-PK	Ø10.0 x 25 mm
	PIS1030-PK	Ø10.0 x 30 mm

Instrument

Item	Item Code	Description
	PIS-0430-DRIVER	Peek interference screw driver



OSTEONIC Co., Ltd.

1206Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea

Tel. +82 2 6902 8400

Fax. +82 2 6902 8401

e-mail. sales@osteonic.com

www.osteonic.com